

HT HT10 HT100

Physical Testing System

SOTAX

Solutions for Pharmaceutical Testing



SOTAX has developed a complete hardness testing line for solid dosage forms. All SOTAX hardness testing systems incorporate the patented Auto-Alignment™ device to ensure proper alignment of any tablet shapes, especially critical shapes such as oblong tablets. All SOTAX hardness testers can be set to test

tablets using "Constant Speed" or "Constant Force" to match a customer's current measurement principles.

▶ The SOTAX HT 1

An ideal solution for IPC, QC and R&D laboratories for testing of single tablets. It can test tablet hardness, thickness and length or diameter as standard with an option for weight measurement.

▶ The SOTAX HT 10



An ideal semi-automated solution for IPC, QC and R&D laboratories for testing of twenty four tablets. It can test tablet hardness, thickness and length or diameter as standard with an option for weight measurement.

▶ The SOTAX HT 100



An ideal fully-automated solution for IPC and QC laboratories where there is a requirement to test large sample volumes. It can test tablet weight, hardness, thickness, width, length or diameter. The HT 100 features a unique system that automatically separates different tablet shapes and sizes without operator intervention - Auto Separator.

CRESCENT
scientific
helping science....



Crescent Scientific Pvt. Ltd.

C-114, Synthofine Industrial Estate, Off Aarey Road,
Goregaon - East, Mumbai - 400 063, Maharashtra, India.

Telephone Numbers +91 - 22 - 2878 7601 / 602

Facsimile Number +91 - 22 - 2875 8852

Email Address crescent@cspl.in web: www.cspl.in