

Encouraged By Your Overwhelming Endorsement
FOR OUR *rapid* RANGE OF EVAPORATORS...



...We Have Expanded Our Footprints

PRESENTING

Microbiology & Laboratory Focused Instruments



FortuneTM

Trusted Quality, Trusted Expertise

FortuneTM Range

Microbiology

Focused Instruments

- Autoclave
- BOD Incubator
- Incubator
- Laminar Airflow
- Oil Bath
- Biosafety Cabinet

Laboratory

Focused Instruments

- Deep Freezer
- Dehumidifier
- Environmental Test Chamber
- Hot Air Oven
- Hot Plate
- Melting cum Boiling Point Apparatus
- Muffle Furnace
- Orbital Shaking Incubators
- Tray Drier
- Ultra-Cryostat
- Water Bath
- Water Bath Incubator Shaker
- Plant Growth Chamber

Trusted Product Quality



Fortune Range is manufactured by a team of professionals with proven experience & expertise in manufacturing microbiology & lab focused instruments since over two decades.

Trusted Service Quality



Fortune Range backed by highly experienced and trained Tech Support Team based at all key locations in India i.e Mumbai, Delhi, Chennai, Bangalore, Hyderabad & Ahmedabad, Providing quick & reliable service.

Trusted Range



Fortune Range provides a **One stop shop** for your all your needs for microbiology & laboratory focused instruments.

Trusted Expertise



Fortune range is backed by high quality experts who can provide consultation related to Microbiology, Lab set-up, Instruments.

BOD Incubator



Application : Microbiology Labs.

Similar Names : Cooling Incubator

Technical Specification / Description :

- Temperature range 50°C to 250 °C.
- Double wall inner Stainless Steel and outer mild steel with two air circulating fans, Microprocessor based Electronic Digital Temperature Controller cum Indicator gives an accuracy of 1°C.
- Acrylic glass powder coat finish door is provided inside.
- Castor wheels are provided at the bottom for easy movement.
- Refrigeration system (CFC free) is provided at the top or bottom

Note : "Rexnord" motors used for internal air circulation.

Optional :

- Touch panel multipurpose temperature controller with recording and PC software
- Other specific size with additional features can be manufactured on request.
- **Validation Port**
Validation at site with desired number of probes for desired duration and modes by the customer.
- IQ/OQ/PQ documents
- Interior illumination and cyclic timer
- GMP Model
- Accuracy of ± 0.1 °C

Technical Specifications

SIZE	4 ft ³	6 ft ³	10 ft ³	12 ft ³
INTERNAL SIZE (H*W*D)	610*455*410 mm	830*505*415 mm	875*570*550 mm	900*650*580 mm
EXTERNAL SIZE (H*W*D)	1090*610*610 mm	1390*660*610 mm	1430*730*760 mm	1460*810*790 mm
CAPACITY (Lts)	120	180	300	360
COMPRESSOR	Yes	Yes	Yes	Yes
COMPRESSOR MOTOR	Yes	Yes	Yes	Yes
CONDENSER	Row	Row	Row	Row
EVAPORATOR	Row	Row	Row	Row
CIRCULATION	2 Nos	2 Nos	2 Nos	2 Nos
CONTROLLER	ON /OFF	ON/OFF	ON/OFF	ON/OFF
SENSOR	PT-100	PT-100	PT-100	PT-100
NO OF TRAYS	2 NOS	3 NOS	3 NOS	4 NOS
SAFETY DOOR	6mm Thick	8mm Thick	8mm Thick	10mm Thick
ILLUMINATION	Yes	Yes	Yes	Yes
HEATER	100 Watts	100 Watts	100 Watts	100 Watts
MAINS	220 VAC	220VAC	220VAC	220VAC

Microbiology & Laboratory Focused Instruments



FortuneTM
Trusted Quality, Trusted Expertise

CONTACT US :

■ Regional Office - Mumbai

A-222, Synthofine Industrial Estate,
Off Aarey Road, Goregaon (East),
Mumbai 400063, Maharashtra, India
Tel: (+91) 22 29275832 / 49726273
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

■ Delhi

713 Vishwadeep Tower,
District Centre, Janakpuri,
New Delhi 110058, India
Tel: (+91) 011 41018861
(+91) 9872927145

Name: Mr. Raghuwinder Yadav
Email: raghu@cspl.in
crescent@cspl.in

■ International Sales

C-114, Synthofine Industrial Estate,
Off Aarey Road, Goregaon East,
Mumbai 400063, Maharashtra, India.
Tel: 29275835 / (+91) 9999162535

Name: Mr. Abhijeet Anand
Email: abhijeet@cspl.in
crescent@cspl.in

■ Chennai

12/202, Ground Floor, 52nd Street,
7th Avenue, Ashok Nagar,
Chennai 600083, India
Tel: (+91) 044 45022263
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Bengaluru

1111/ 1112, 'B' Wing, 11th Floor,
Mittal Towers, No.6, M G Road,
Bengaluru 560001, Karnataka, India
Tel: (+91) 080 41120733
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Hyderabad

201, Lingapur House,
Amrutha Estates, Himayatnagar,
Hyderabad 500 029,
Telangana, India
Tel: (+91) 040 66848836 / 40102005
(+91) 9949651182

Name: Mr. Ganesh .M
Email: ganesh@cspl.in
crescent@cspl.in

■ Ahmedabad

709, 7th Floor, Shitiratna Building,
Opp. Panchwati Circle, Ambawadi,
Ellis Bridge, Ahmedabad - 380006
Tel: 079-26444180 / 40029180
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

www.cspl.in