

Encouraged By Your Overwhelming Endorsement
FOR OUR *rapid* RANGE OF EVAPORATORS...



...We Have Expanded Our Footprints

PRESENTING

Microbiology & Laboratory Focused Instruments



FortuneTM

Trusted Quality, Trusted Expertise

FortuneTM Range

Microbiology

Focused Instruments

- Autoclave
- BOD Incubator
- Incubator
- Laminar Airflow
- Oil Bath
- Biosafety Cabinet

Laboratory

Focused Instruments

- Deep Freezer
- Dehumidifier
- Environmental Test Chamber
- Hot Air Oven
- Hot Plate
- Melting cum Boiling Point Apparatus
- Muffle Furnace
- Orbital Shaking Incubators
- Tray Drier
- Ultra-Cryostat
- Water Bath
- Water Bath Incubator Shaker
- Plant Growth Chamber

Trusted Product Quality



Fortune Range is manufactured by a team of professionals with proven experience & expertise in manufacturing microbiology & lab focused instruments since over two decades.

Trusted Service Quality



Fortune Range backed by highly experienced and trained Tech Support Team based at all key locations in India i.e Mumbai, Delhi, Chennai, Bangalore, Hyderabad & Ahmedabad, Providing quick & reliable service.

Trusted Range



Fortune Range provides a **One stop shop** for your all your needs for microbiology & laboratory focused instruments.

Trusted Expertise



Fortune range is backed by high quality experts who can provide consultation related to Microbiology, Lab set-up, Instruments.

Deep Freezer

-20°C Deep Freezer

VERTICAL MODEL



Application : Used in all labs for storing samples, Medicines, Injectables, Food samples & Blood samples etc.

Technical Specification / Description :

- These deep freezers are designed for quick-freezing and storing of blood ITS components, serum, vaccines, biological and medical specimens, clinical samples at low temperature.
- These freezers are with specially engineered superior insulation and specially designed refrigeration system constructed of double walls.
- The exterior is made of sheet steel while the interior of Stainless Steel having puff insulation to minimize temperature loss.
- Capillary controlled balanced hermetically sealed refrigeration is most economical, trouble free and exceptionally quiet.

HORIZONTAL MODEL



Note :
Automatic Voltage Stabilizer with time delay rating : 2 KVA for the above to be provided by the customer.

Optional :

- 40°C deep freezer is also available.

Technical Specifications

	170 LTS	265 LTS
CAPACITY	170 LTS	265 LTS
EXTERNAL SIZE (W*D*H)	80 X 60 X 147 cm	115 X 60 X 147 cm
INTERNAL SIZE (W*D*H)	60 X 40 X 70 cm	95 X 40 X 70 cm
COMPRESSOR	1.5 TON	2 TON
COIL	Copper Coil	Copper Coil
COMPRESSOR MOTOR	1/5 HP	1/5 HP
FAN BLADE	LEAF	LEAF
INSULATION	PUFF	PUFF
CONTROLLER	ON / OFF	ON / OFF
SENSOR	PT -100	PT -100
SWITCH	10 AMPS ON / OFF	10 AMPS ON / OFF
TOTAL POWER	15 AMPS	15 AMPS
MAINS	220V Single Phase	220V Single Phase

Microbiology & Laboratory Focused Instruments



FortuneTM
Trusted Quality, Trusted Expertise

CONTACT US :

■ Regional Office - Mumbai

A-222, Synthofine Industrial Estate,
Off Aarey Road, Goregaon (East),
Mumbai 400063, Maharashtra, India
Tel: (+91) 22 29275832 / 49726273
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

■ Delhi

713 Vishwadeep Tower,
District Centre, Janakpuri,
New Delhi 110058, India
Tel: (+91) 011 41018861
(+91) 9872927145

Name: Mr. Raghuwinder Yadav
Email: raghu@cspl.in
crescent@cspl.in

■ International Sales

C-114, Synthofine Industrial Estate,
Off Aarey Road, Goregaon East,
Mumbai 400063, Maharashtra, India.
Tel: 29275835 / (+91) 9999162535

Name: Mr. Abhijeet Anand
Email: abhijeet@cspl.in
crescent@cspl.in

■ Chennai

12/202, Ground Floor, 52nd Street,
7th Avenue, Ashok Nagar,
Chennai 600083, India
Tel: (+91) 044 45022263
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Bengaluru

1111/ 1112, 'B' Wing, 11th Floor,
Mittal Towers, No.6, M G Road,
Bengaluru 560001, Karnataka, India
Tel: (+91) 080 41120733
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Hyderabad

201, Lingapur House,
Amrutha Estates, Himayatnagar,
Hyderabad 500 029,
Telangana, India
Tel: (+91) 040 66848836 / 40102005
(+91) 9949651182

Name: Mr. Ganesh .M
Email: ganesh@cspl.in
crescent@cspl.in

■ Ahmedabad

709, 7th Floor, Shitiratna Building,
Opp. Panchwati Circle, Ambawadi,
Ellis Bridge, Ahmedabad - 380006
Tel: 079-26444180 / 40029180
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

www.cspl.in