

Encouraged By Your Overwhelming Endorsement  
**FOR OUR *rapid* RANGE OF EVAPORATORS...**



**...We Have Expanded Our Footprints**

# PRESENTING

## Microbiology & Laboratory Focused Instruments



# Fortune<sup>TM</sup>

Trusted Quality, Trusted Expertise

## Fortune<sup>TM</sup> Range

### Microbiology

#### Focused Instruments

- Autoclave
- BOD Incubator
- Incubator
- Laminar Airflow
- Oil Bath
- Biosafety Cabinet

### Laboratory

#### Focused Instruments

- Deep Freezer
- Dehumidifier
- Environmental Test Chamber
- Hot Air Oven
- Hot Plate
- Melting cum Boiling Point Apparatus
- Muffle Furnace
- Orbital Shaking Incubators
- Tray Drier
- Ultra-Cryostat
- Water Bath
- Water Bath Incubator Shaker
- Plant Growth Chamber

## Trusted Product Quality



Fortune Range is manufactured by a team of professionals with proven experience & expertise in manufacturing microbiology & lab focused instruments since over two decades.

## Trusted Service Quality



Fortune Range backed by highly experienced and trained Tech Support Team based at all key locations in India i.e Mumbai, Delhi, Chennai, Bangalore, Hyderabad & Ahmedabad, Providing quick & reliable service.

## Trusted Range



Fortune Range provides a **One stop shop** for your all your needs for microbiology & laboratory focused instruments.

## Trusted Expertise



Fortune range is backed by high quality experts who can provide consultation related to Microbiology, Lab set-up, Instruments.

## Hot Air Oven



**Application :** Used to dry heat and sterilize the Lab products in Microbiology Labs, in Pharma, Food, Chemical and many other Industries.

**Similar Names :** Drier

### Technical Specification / Description :

- Temperature Range 180°C. Digital Temperature Controller with an Accuracy of  $\pm 1^\circ\text{C}$ .
- Heating Element are provided at the bottom only.
- Double walled inner stainless steel and outside mild steel powder coat finish with trays and with ,air circulating fan works on 220 Volts AC.
- LCD/Touch panel micro PLC controller with recordings, software and Printer

### GMP Model

- Elements are placed at the sides and bottom.
- Triple walled inner and outer Stainless Steel 304 grade finish with shelves and air circulating fan works on 220 Volts AC.
- Digital Temperature Controller Micro Processor Based is Provided with an accuracy of  $\pm 2^\circ\text{C}$ .
- Temperature  $5^\circ\text{C}$  above ambient to  $220^\circ\text{C}$  but recommended working temperature is  $180^\circ\text{C}$ .

### Technical Specifications

<b>INTERNAL SIZE (W*D*H)</b>	355 X 355 X 355 mm	455 X 455D X 455 mm	455 X 455 X 605 mm	605 X 605 X 605 mm	605 X 605 X 910 mm
<b>EXTERNAL SIZE (W*D*H)</b>	630 X 560 X 755 mm	730 X 660 X 960 mm	730 X 660 X 1080 mm	880 X 810 X 1080 mm	880 X 810 X 1390 mm
<b>CAPACITY</b>	1.6 ft <sup>3</sup>	3.4 ft <sup>3</sup>	4.8 ft <sup>3</sup>	7.9 ft <sup>3</sup>	11.9 ft <sup>3</sup>
<b>HEATER RATING</b>	1500 W	1800 W	2000 W	2400 W	3000 W
<b>MOTOR</b>	1/50 HP	1/50 HP	1/50 HP	1/50 HP	1/12 HP
<b>BLOWER</b>	5*2 ACW	5*2 ACW	5*2 ACW	6*2 ACW	7*2 CW
<b>NO OF TRAYS</b>	2	2	3	3	4
<b>DOOR LOCK</b>	DSNL	DSNL	DSNL	DSNL	DSNL
<b>ON / OFF SWITCH</b>	12 AMPS	12 AMPS	12 AMPS	12 AMPS	16 AMPS
<b>ROCKER SWITCH</b>	12 AMPS	12 AMPS	12 AMPS	12 AMPS	12 AMPS
<b>CONTROLLER</b>	PID (ON-OFF)	PID(ON-OFF)	PID (ON-OFF)	PID (ON-OFF)	PID (ON-OFF)
<b>POWER CORD</b>	6 AMPS	15 AMPS	15 AMPS	15 AMPS	20 AMPS

## Microbiology & Laboratory Focused Instruments



**Fortune**<sup>TM</sup>  
Trusted Quality, Trusted Expertise

### CONTACT US :

#### ■ Regional Office - Mumbai

A-222, Synthofine Industrial Estate,  
Off Aarey Road, Goregaon (East),  
Mumbai 400063, Maharashtra, India  
Tel: (+91) 22 29275832 / 49726273  
(+91) 9867599915

Name: Mr. Hiten Rana  
Email: [hiten@cspl.in](mailto:hiten@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Delhi

713 Vishwadeep Tower,  
District Centre, Janakpuri,  
New Delhi 110058, India  
Tel: (+91) 011 41018861  
(+91) 9872927145

Name: Mr. Raghuwinder Yadav  
Email: [raghu@cspl.in](mailto:raghu@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ International Sales

C-114, Synthofine Industrial Estate,  
Off Aarey Road, Goregaon East,  
Mumbai 400063, Maharashtra, India.  
Tel: 29275835 / (+91) 9999162535

Name: Mr. Abhijeet Anand  
Email: [abhijeet@cspl.in](mailto:abhijeet@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Chennai

12/202, Ground Floor, 52nd Street,  
7th Avenue, Ashok Nagar,  
Chennai 600083, India  
Tel: (+91) 044 45022263  
(+91) 9884266900

Name: Mr. G. Palaniappan  
Email: [palani@cspl.in](mailto:palani@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Bengaluru

# 1111/ 1112, 'B' Wing, 11th Floor,  
Mittal Towers, No.6, M G Road,  
Bengaluru 560001, Karnataka, India  
Tel: (+91) 080 41120733  
(+91) 9884266900

Name: Mr. G. Palaniappan  
Email: [palani@cspl.in](mailto:palani@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Hyderabad

201, Lingapur House,  
Amrutha Estates, Himayatnagar,  
Hyderabad 500 029,  
Telangana, India  
Tel: (+91) 040 66848836 / 40102005  
(+91) 9949651182

Name: Mr. Ganesh .M  
Email: [ganesh@cspl.in](mailto:ganesh@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Ahmedabad

709, 7th Floor, Shitiratna Building,  
Opp. Panchwati Circle, Ambawadi,  
Ellis Bridge, Ahmedabad - 380006  
Tel: 079-26444180 / 40029180  
(+91) 9867599915

Name: Mr. Hiten Rana  
Email: [hiten@cspl.in](mailto:hiten@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

[www.cspl.in](http://www.cspl.in)