

Encouraged By Your Overwhelming Endorsement  
**FOR OUR *rapid* RANGE OF EVAPORATORS...**



**...We Have Expanded Our Footprints**

# PRESENTING

## Microbiology & Laboratory Focused Instruments



# Fortune<sup>TM</sup>

Trusted Quality, Trusted Expertise

## Fortune<sup>TM</sup> Range

### Microbiology

#### Focused Instruments

- Autoclave
- BOD Incubator
- Incubator
- Laminar Airflow
- Oil Bath
- Biosafety Cabinet

### Laboratory

#### Focused Instruments

- Deep Freezer
- Dehumidifier
- Environmental Test Chamber
- Hot Air Oven
- Hot Plate
- Melting cum Boiling Point Apparatus
- Muffle Furnace
- Orbital Shaking Incubators
- Tray Drier
- Ultra-Cryostat
- Water Bath
- Water Bath Incubator Shaker
- Plant Growth Chamber

## Trusted Product Quality



Fortune Range is manufactured by a team of professionals with proven experience & expertise in manufacturing microbiology & lab focused instruments since over two decades.

## Trusted Service Quality



Fortune Range backed by highly experienced and trained Tech Support Team based at all key locations in India i.e Mumbai, Delhi, Chennai, Bangalore, Hyderabad & Ahmedabad, Providing quick & reliable service.

## Trusted Range



Fortune Range provides a **One stop shop** for your all your needs for microbiology & laboratory focused instruments.

## Trusted Expertise



Fortune range is backed by high quality experts who can provide consultation related to Microbiology, Lab set-up, Instruments.

## Tray Drier-Double Side Door

**Application :** Agricultural Drying

**Similar Names :** Industrial Drier

### Technical Specification / Description :

- The drier is double walled with insulation and is fitted with best quality heaters.
- A suitable blower driven by motor circulates the hot air.
- The temperature is controlled by a Micro Processor based Digital Temperature Controller cum Indicator.
- The unit is supplied complete with panel box having indication lamps, temperature indicator, starter, contactor etc. The unit is supplied without trolley and with trays.
- SS Interior and MS Exterior nicely painted.



<b>NO OF TRAYS</b>	12 Trays	24 Trays	48 Trays
<b>MOTOR</b>	1/4 HP, 3 Phase	1/2 HP, 3 Phase	1 HP Motor, 3 Phase
<b>HEATING LOAD</b>	3 KW	6 KW	9 KW
<b>TEMPERATURE</b>	100°C	100°C	100°C

### ADDITIONAL TROLLEY AND TRAYS AT EXTRA COST

- FOR SS / ALUMINIUM TRAYS 20 SWG  
SIZE: 32 X 16 X 1.25"
- DIGITAL TEMPERATURE CONTROLLER  
CUM INDICATOR
- FOR TROLLEY MILD STEEL / SS

## Microbiology & Laboratory Focused Instruments



**Fortune**<sup>TM</sup>  
Trusted Quality, Trusted Expertise

### CONTACT US :

#### ■ Regional Office - Mumbai

A-222, Synthofine Industrial Estate,  
Off Aarey Road, Goregaon (East),  
Mumbai 400063, Maharashtra, India  
Tel: (+91) 22 29275832 / 49726273  
(+91) 9867599915

Name: Mr. Hiten Rana  
Email: [hiten@cspl.in](mailto:hiten@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Delhi

713 Vishwadeep Tower,  
District Centre, Janakpuri,  
New Delhi 110058, India  
Tel: (+91) 011 41018861  
(+91) 9872927145

Name: Mr. Raghuwinder Yadav  
Email: [raghu@cspl.in](mailto:raghu@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ International Sales

C-114, Synthofine Industrial Estate,  
Off Aarey Road, Goregaon East,  
Mumbai 400063, Maharashtra, India.  
Tel: 29275835 / (+91) 9999162535

Name: Mr. Abhijeet Anand  
Email: [abhijeet@cspl.in](mailto:abhijeet@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Chennai

12/202, Ground Floor, 52nd Street,  
7th Avenue, Ashok Nagar,  
Chennai 600083, India  
Tel: (+91) 044 45022263  
(+91) 9884266900

Name: Mr. G. Palaniappan  
Email: [palani@cspl.in](mailto:palani@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Bengaluru

# 1111/ 1112, 'B' Wing, 11th Floor,  
Mittal Towers, No.6, M G Road,  
Bengaluru 560001, Karnataka, India  
Tel: (+91) 080 41120733  
(+91) 9884266900

Name: Mr. G. Palaniappan  
Email: [palani@cspl.in](mailto:palani@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Hyderabad

201, Lingapur House,  
Amrutha Estates, Himayatnagar,  
Hyderabad 500 029,  
Telangana, India  
Tel: (+91) 040 66848836 / 40102005  
(+91) 9949651182

Name: Mr. Ganesh .M  
Email: [ganesh@cspl.in](mailto:ganesh@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

#### ■ Ahmedabad

709, 7th Floor, Shitiratna Building,  
Opp. Panchwati Circle, Ambawadi,  
Ellis Bridge, Ahmedabad - 380006  
Tel: 079-26444180 / 40029180  
(+91) 9867599915

Name: Mr. Hiten Rana  
Email: [hiten@cspl.in](mailto:hiten@cspl.in)  
[crescent@cspl.in](mailto:crescent@cspl.in)

[www.cspl.in](http://www.cspl.in)