

Encouraged By Your Overwhelming Endorsement
FOR OUR *rapid* RANGE OF EVAPORATORS...



...We Have Expanded Our Footprints

PRESENTING

Microbiology & Laboratory Focused Instruments



FortuneTM

Trusted Quality, Trusted Expertise

FortuneTM Range

Microbiology

Focused Instruments

- Autoclave
- BOD Incubator
- Incubator
- Laminar Airflow
- Oil Bath
- Biosafety Cabinet

Laboratory

Focused Instruments

- Deep Freezer
- Dehumidifier
- Environmental Test Chamber
- Hot Air Oven
- Hot Plate
- Melting cum Boiling Point Apparatus
- Muffle Furnace
- Orbital Shaking Incubators
- Tray Drier
- Ultra-Cryostat
- Water Bath
- Water Bath Incubator Shaker
- Plant Growth Chamber

Trusted Product Quality



Fortune Range is manufactured by a team of professionals with proven experience & expertise in manufacturing microbiology & lab focused instruments since over two decades.

Trusted Service Quality



Fortune Range backed by highly experienced and trained Tech Support Team based at all key locations in India i.e Mumbai, Delhi, Chennai, Bangalore, Hyderabad & Ahmedabad, Providing quick & reliable service.

Trusted Range



Fortune Range provides a **One stop shop** for your all your needs for microbiology & laboratory focused instruments.

Trusted Expertise



Fortune range is backed by high quality experts who can provide consultation related to Microbiology, Lab set-up, Instruments.

Ultra Cryostat Circulator



Application : Pharma, Environmental, Chemical and Universities

Similar Names : Chiller

Technical Specification / Description :

- The inner tank shall be made out of Stainless Steel and is insulated with thermorix sheet.
- Unit is mounted on a sturdy iron frame.
- The refrigeration system consists of a compressor, condenser and a fan which are fixed at the bottom of the frame.
- The evaporator coil is placed inside the SS Tank. (R-134A) gas is charged to the refrigeration system which brings down the temperature.
- A stirrer keeps the temperature uniform and accuracy is controlled by a digital Controller with 3.5" digit led with resolution of 0.5°C.
- The complete unit is covered with mild steel sheet and powder coated in a pleasing colour.
- Four wheels are provided for easy movement.

Optional :

- Special Type Magnetic Pumping unit with an output of 10 ltr/min at 2 mtr head to enable water circulation.

- Voltage stabilizer 3 KVA is a must for the above unit to be provided by the customer.

Technical Specifications

TANK CAPACITY	15 Lts
INNER SIZE (W*D*H)	25 X 25 X 30 cm
EXTERNAL SIZE (W*D*H)	380 X 560 X 660 mm
ACCURACY	0.5°
MEDIUM	Water
CONTROLLER	OFF Type
SENSOR	PT -100
CIRCULATION	Magnetic Pump
INSULATION	Thermorexin
POWER CONSUMPTION	10 AMPS

Microbiology & Laboratory Focused Instruments



FortuneTM
Trusted Quality, Trusted Expertise

CONTACT US :

■ Regional Office - Mumbai

A-222, Synthofine Industrial Estate,
Off Aarey Road, Goregaon (East),
Mumbai 400063, Maharashtra, India
Tel: (+91) 22 29275832 / 49726273
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

■ Delhi

713 Vishwadeep Tower,
District Centre, Janakpuri,
New Delhi 110058, India
Tel: (+91) 011 41018861
(+91) 9872927145

Name: Mr. Raghuwinder Yadav
Email: raghu@cspl.in
crescent@cspl.in

■ International Sales

C-114, Synthofine Industrial Estate,
Off Aarey Road, Goregaon East,
Mumbai 400063, Maharashtra, India.
Tel: 29275835 / (+91) 9999162535

Name: Mr. Abhijeet Anand
Email: abhijeet@cspl.in
crescent@cspl.in

■ Chennai

12/202, Ground Floor, 52nd Street,
7th Avenue, Ashok Nagar,
Chennai 600083, India
Tel: (+91) 044 45022263
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Bengaluru

1111/ 1112, 'B' Wing, 11th Floor,
Mittal Towers, No.6, M G Road,
Bengaluru 560001, Karnataka, India
Tel: (+91) 080 41120733
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Hyderabad

201, Lingapur House,
Amrutha Estates, Himayatnagar,
Hyderabad 500 029,
Telangana, India
Tel: (+91) 040 66848836 / 40102005
(+91) 9949651182

Name: Mr. Ganesh .M
Email: ganesh@cspl.in
crescent@cspl.in

■ Ahmedabad

709, 7th Floor, Shitiratna Building,
Opp. Panchwati Circle, Ambawadi,
Ellis Bridge, Ahmedabad - 380006
Tel: 079-26444180 / 40029180
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

www.cspl.in