

Encouraged By Your Overwhelming Endorsement
FOR OUR *rapid* RANGE OF EVAPORATORS...



...We Have Expanded Our Footprints

PRESENTING

Microbiology & Laboratory Focused Instruments



FortuneTM

Trusted Quality, Trusted Expertise

FortuneTM Range

Microbiology

Focused Instruments

- Autoclave
- BOD Incubator
- Incubator
- Laminar Airflow
- Oil Bath
- Biosafety Cabinet

Laboratory

Focused Instruments

- Deep Freezer
- Dehumidifier
- Environmental Test Chamber
- Hot Air Oven
- Hot Plate
- Melting cum Boiling Point Apparatus
- Muffle Furnace
- Orbital Shaking Incubators
- Tray Drier
- Ultra-Cryostat
- Water Bath
- Water Bath Incubator Shaker
- Plant Growth Chamber

Trusted Product Quality



Fortune Range is manufactured by a team of professionals with proven experience & expertise in manufacturing microbiology & lab focused instruments since over two decades.

Trusted Service Quality



Fortune Range backed by highly experienced and trained Tech Support Team based at all key locations in India i.e Mumbai, Delhi, Chennai, Bangalore, Hyderabad & Ahmedabad, Providing quick & reliable service.

Trusted Range



Fortune Range provides a **One stop shop** for your all your needs for microbiology & laboratory focused instruments.

Trusted Expertise



Fortune range is backed by high quality experts who can provide consultation related to Microbiology, Lab set-up, Instruments.

Water Bath Incubator Shaker

Application : Biotechnology, Clinical, Micro Biology, Medical, Environmental or Industrial Applications.

Similar Names : Metabolic Shaker

Technical Specification / Description :

- Temperature range from room temperature to 100°C.
- Digital temperature control with an accuracy of $\pm 1^\circ\text{C}$.
- Double walled inside made of Stainless Steel 304 quality and outside mild steel sheet powder coat painted.
- Oscillating tray is riding on ball rollers.
- Shaking speed range 40 to 140 rpm. Built-in water drain and digital rpm indicator provided.
- Works on 220/230 volts AC.



SHAKING RACKS MADE OF STAINLESS STEEL FOR TYPE ISB 1 & ISB 2 - To Hold

16 NOS.	25 ML	Conical Flask
2 NOS.	50 ML	Conical Flask
9 NOS.	100 ML	Conical Flask
5 NOS.	250 ML	Conical Flask
4 NOS.	500 ML	Conical Flask
2 NOS.	1000 ML	Conical Flask
25 NOS.	25 MM	Dia Test Tube

SHAKING RACKS MADE OF STAINLESS STEEL FOR TYPE ISB-3 - To Hold

20 NOS.	25 ML	Conical Flask
16 NOS.	50 ML	Conical Flask
12 NOS.	100 ML	Conical Flask
9 NOS.	250 ML	Conical Flask
6 NOS.	500 ML	Conical Flask
4 NOS.	1000 ML	Conical Flask
35 NOS.	25 MM	Dia Test Tube

SHAKING TRAY SIZE: INSIDE (STAINLESS STEEL 304 QTY):

ISB-1	275 X 275 X 100 MM
ISB-2	275 X 275 X 170 MM
ISB-3	405 X 300 X 170 MM

Optional :

- GMP model
- IQ/OQ/PQ documents.

Technical Specifications

INTERNAL SIZE	27.5" X 27.5" X 7.5"	27.5" X 27.5" X 15"	40.5" X 30" X 15"
	275 X 275 X 75 (mm)	275 X 275 X 150 (mm)	405 X 300 X 150 (mm)
EXTERNAL SIZE (W*H*D)	380*580*200 mm	380*580*230 mm	510*600*230 mm
TEMPERATURE RANGE	Ambient To 100 ^o	Ambient To 100 ^o	Ambient To 100 ^o
HEATER	1500 WATTS	2500 WATTS	2500 WATTS
CONTROLLER	Micro Processor PID	Micro Processor PID	Micro Processor PID
SENSOR	PT-100	PT-100	PT-100
RPM RANGE	0-140 RPM	0-140 RPM	0-140 RPM
INSULATION	Mineral Wool	Mineral Wool	Mineral Wool
SPEED CONTROL	By Regulator	By Regulator	By Regulator
MAINS	220 VAC Single Phase	220 VAC Single Phase	220 VAC Single Phase

Microbiology & Laboratory Focused Instruments



FortuneTM
Trusted Quality, Trusted Expertise

CONTACT US :

■ Regional Office - Mumbai

A-222, Synthofine Industrial Estate,
Off Aarey Road, Goregaon (East),
Mumbai 400063, Maharashtra, India
Tel: (+91) 22 29275832 / 49726273
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

■ Delhi

713 Vishwadeep Tower,
District Centre, Janakpuri,
New Delhi 110058, India
Tel: (+91) 011 41018861
(+91) 9872927145

Name: Mr. Raghuwinder Yadav
Email: raghu@cspl.in
crescent@cspl.in

■ International Sales

C-114, Synthofine Industrial Estate,
Off Aarey Road, Goregaon East,
Mumbai 400063, Maharashtra, India.
Tel: 29275835 / (+91) 9999162535

Name: Mr. Abhijeet Anand
Email: abhijeet@cspl.in
crescent@cspl.in

■ Chennai

12/202, Ground Floor, 52nd Street,
7th Avenue, Ashok Nagar,
Chennai 600083, India
Tel: (+91) 044 45022263
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Bengaluru

1111/ 1112, 'B' Wing, 11th Floor,
Mittal Towers, No.6, M G Road,
Bengaluru 560001, Karnataka, India
Tel: (+91) 080 41120733
(+91) 9884266900

Name: Mr. G. Palaniappan
Email: palani@cspl.in
crescent@cspl.in

■ Hyderabad

201, Lingapur House,
Amrutha Estates, Himayatnagar,
Hyderabad 500 029,
Telangana, India
Tel: (+91) 040 66848836 / 40102005
(+91) 9949651182

Name: Mr. Ganesh .M
Email: ganesh@cspl.in
crescent@cspl.in

■ Ahmedabad

709, 7th Floor, Shitiratna Building,
Opp. Panchwati Circle, Ambawadi,
Ellis Bridge, Ahmedabad - 380006
Tel: 079-26444180 / 40029180
(+91) 9867599915

Name: Mr. Hiten Rana
Email: hiten@cspl.in
crescent@cspl.in

www.cspl.in